

RNA-Seq Analysis Troubleshooting

RNA-seq Workflow







All NGS analyses require that the quality of the raw data is assessed prior to any downstream analysis.

The quality checks at this stage in the workflow include:

- 1. Checking the **quality of the base calls** to ensure that there were no issues during sequencing
- 2. Examining the reads to ensure their quality metrics adhere to our expectations for our experiment
- 3. Exploring reads for contamination





Sequence	Count	Percentage	Possible Sou
000000000000000000000000000000000000000	1838931	25.27061807325126	No Hit
ADDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	22246	0.3057048740042707	No Hit
одерозоророзорозорозорозорозорозорозорозор	19143	0.26306340029954844	No Hit
204202020202020202020202020202020202020	14083	0.1935288025084125	No Hit
GGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGGG	13303	0.18281003051689354	No Hit
AKAKAKDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	12912	0.17743690250576033	No Hit
OCCOCCOCOCOCOCOCOCOCOCOCOCOCOCOCOCOCOCOC	11561	0.1588714397358345	No Hit
000000000000000000000000000000000000000	11428	0.15704375169112678	No Hit
AAAAAAAAAAAAADDDDDDDDDDDDDDDDDDDDDDDDDD	11152	0.1532509554479739	No Hit
GCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC	9922	0.1363482765382709	No Hit
${\color{blue}AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA$	9693	0.1332013550176839	No Hit
AKDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	9340	0.12835042359075288	No Hit
AKKKKASDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	8557	0.11759042555311268	No Hit
AAAAADDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	8156	0.11207987738824203	No Hit
DDDDDADDDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	7471	0.10266659685722856	No Hit
DDDDDDADDDDDDDDDDDDDDDDDDDDDDDDDDDDDDD	7294	0.10023426013607618	No Hit

Troubleshooting raw data quality problems:

- Poor quality data (problems at sequencing facility)
 - Poor quality across sequence
 - Drop in quality in the middle
 - Large percentage of sequences with low mean quality scores

Troubleshooting raw data quality problems:

- Issues based on read sequence expectations
 - Unexpected %GC for organism or % of each nucleotide does not remain similar across the read (except for first 10-12 bases)
 - Contaminating sequences: different species, adapters, vector, mitochondrial/rRNA
 - High level of sequence duplications
 - low complexity library, too many cycles of PCR amplification / too little starting material
 - Over-represented sequences more than 1-2%, unless expected based on experimental design
 - contaminating sequences: adapters, vector, mitochondrial/rRNA


Raw Data QC Goals:

- Identify sequencing problems and determine whether there is a need to contact the sequencing facility
- Identify over-represented contaminating sequences
- Gain insight into library complexity (rRNA contamination, duplications)
- Ensure organism is properly represented by %GC content

Raw Data QC Goals:

Can we identify a degraded RNA-Seq sample (low RIN #) using these raw data QC metrics?

Since reads from degraded samples are generally just shorter, the quality of the sequenced nucleotides should be fine. At this step, degraded libraries will not likely affect the quality metrics.



Evaluating the **quality of the aligned data** can give important information about the quality of the library. The quality checks at this stage in the workflow include:

- 1. Checking the total percent of reads aligning to the genome
- 2. Determining the percent uniquely mapping reads
- 3. Examining the total number of reads aligning to each sample
- 4. Checking percent of paired-end reads that are properly paired

Troubleshooting aligned data quality problems:

- Low percentage (< 80%) of reads aligned
 - poor quality reads, contaminating sequences, inappropriate alignment parameters chosen, inappropriate reference genome chosen, poor quality reference genome
- Low percentage (< 60%) of uniquely aligning reads
 - low number of total reads aligning, organism has high number of paralogous genes, very short read length, low quality bases
- Large differences in sequencing depth between samples
 - library prep / sequencing
- For paired-end data: large number of reads not properly paired
 - poor quality reads

Aligned Data QC Goals:

- Ensure the library depth and percentage of reads mapping to each sample is similar
- Identify poor alignment parameters or low quality library
- Discover contamination from another organism

Quality Checks: Quantified Data

